




# *CHALLENGE 4.*

# *TEMPORAL SIGNALS*


Debora Gil

# Starting Point





Detection of epileptic seizures in EEG signals


# Seizure Patterns in EEG




normal




seizure



# System Pipeline





**Data Preparation.** Extract temporal windows from EEG signals



**Window Seizure Detection.** Identify seizure presence in Windows from the analysis of NCh EEG signals


# Window Seizure Detection



## Classification Approach




## Anomaly Detector




# AI Challenges

## 1. How to combine EEG electrode signals (Channel Fusion)





## 2. How to include temporal information of EEG signals





# Channel Fusion


Fusion of EEG channels at input level (like image CNNs)




Fusion of EEG channels at feature level (after extracting relevant features for each sensor)

# Temporal Information

## 2. How to include temporal information of EEG signals



Recurrent Networks (LSTM)




Transformers



# Temporal Information


Combine channels using temporal information



$X=[21,T]$  is considered a temporal signal of length  $T=128$  represented by the 21 channels

Process consecutive temporal windows to classify the last one

BackBone network to extract LSTM input features



# AI Challenges

## 3. Personalization level of models

Population Approach.




Test: unseen subjects

Personalized Approach.



Test: unseen event of same set of subjects

# Challenge PipeLine



# *Specific Objectives*

---

## 1. Implement two systems:

1.1 Baseline single window classifier with Channel fusion

1.2 Temporal approach using LSTM

## 2. Compare performance at different generalization levels:

2.1 Personalized Model

2.2 Poblational Model

# Session Schedule

desembre	dimarts		
	2	Challenge 4: Detection of Epileptic Seizure in EEG. Challenge presentation. Data Base Exploring Materials in VC: Challenge Presentation. CHB-MIT Dataset.	
	divendres		
	5	Base Line Architecture for Channel Fusion. Materials in VC: Architectures.Channel Fusion	
	dimarts		
	9	Base Line Architecture for Channel Fusion. Follow-up Materials in VC: Architectures.Channel Fusion	
	divendres		
	12	Experimental Design. Train and Test of Base Line Architecture at two different generalization levels. Materials in VC: Validation&Verification	
	dimarts		
	16	Recurrent Architecture for Processing of Temporal Information. Materials in VC: Architectures. Processing Temporal Data	
	divendres		
	19	Recurrent Architecture for Processing of Temporal Information. Follow-up Materials in VC: Architectures.Processing Temporal Data	
			Gener-Febrer
			Defense of the Challenge (Check Calendar of Exams on the Web Site)

---



*Dataset*



# CHB-MIT DataSet



Recordings at 256 MHz from 24 subjects with several seizures. Seizure start-end time is annotated.

EEG device with 21 electrodes



# Annotated windows

Signals resampled at 128 MHz and cropped in 1sec temporal windows





Windows annotations:  
0-normal  
1-ictal

pre and post ictal intervals have to be discarded

# Annotated windows

Ictal windows with 80% overlap for ictal data augmentation



# Annotated windows files

---

## 1. npz files with EEG Windows signals:

PatID\_seizure\_EEGwindow\_1.npz

PatID = {chb\_01, ....., chb24}

EEG\_win is ndarray of size [NWinPatID, 21, 128]

Number of annotated windows



The diagram consists of three arrows pointing upwards towards the dimensions of the ndarray [NWinPatID, 21, 128]. The first arrow originates from the text 'Number of annotated windows' and points to 'NWinPatID'. The second arrow originates from the text 'Number of EEG Channels' and points to '21'. The third arrow originates from the text 'Length of the 1 sec. window' and points to '128'.

Number of EEG Channels

Length of the 1 sec. window


# Annotated windows files

## 2. parquet with EEG Windows Metadata: PatID\_seizure\_metadata\_1.parquet

Dataframe with class, recording identifier (.edf), normal/ictal Interval identifiers of each window

	class	filename_interval	global_interval	filename
0	0	1	1	chb01_03.edf
1	0	1	1	chb01_03.edf
2	0	1	1	chb01_03.edf
...	...	...	...	...
26525	1	2	7	chb01_26.edf
26526	1	2	7	chb01_26.edf
26527	1	2	7	chb01_26.edf
26528	1	2	7	chb01_26.edf

Usefull for data split  
at the two levels of  
personalization




```
metadata=pd.read_parquet(
file,engine='fastparquet')
```

# Annotated windows metadata


Personalized model  
using seizure Interval  
split

chb01\_03.edf


normal1



seizure1




normal2




***global\_Interval*** identifier for all recordings of a given patient

chb01\_26.edf


normal21




seizure6



normal22



seizure7






# Annotated windows metadata


Personalized model  
using patient recording  
split

chb01\_03.edf


normal1



seizure1




normal2




***filename*** identifies different recordings of the patient

chb01\_26.edf


normal21




seizure6



normal22



seizure7



# Annotated windows metadata

Population model using patient split  
(leave-one-out)

***filename.Split('\_')[0]***  
identifies the patient

	class	filename_interval	global_interval	filename
0	0	1	1	chb01_03.edf
1	0	1	1	chb01_03.edf
2	0	1	1	chb01_03.edf
...	...	...	...	
26525	1	2	7	chb01_26.edf
26526	1	2	7	chb01_26.edf
26527	1	2	7	chb01_26.edf
26528	1	2	7	chb01_26.edf

# *CHALLENGE 4.*

# *TEMPORAL SIGNALS*

Debora Gil